



Cisco HealthPresence Solution Design Guide Addendum

Cisco HealthPresence Version 2.0 Addendum for Neurosynaptic Telemedicine
Medical Devices
November 18, 2011

Americas Headquarters
Cisco Systems, Inc.
170 West Tasman Drive
San Jose, CA 95134-1706
USA
<http://www.cisco.com>
Tel: 408 526-4000
800 553-NETS (6387)
Fax: 408 527-0883

Text Part Number: OL-25211-01

THE SPECIFICATIONS AND INFORMATION REGARDING THE PRODUCTS IN THIS MANUAL ARE SUBJECT TO CHANGE WITHOUT NOTICE. ALL STATEMENTS, INFORMATION, AND RECOMMENDATIONS IN THIS MANUAL ARE BELIEVED TO BE ACCURATE BUT ARE PRESENTED WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED. USERS MUST TAKE FULL RESPONSIBILITY FOR THEIR APPLICATION OF ANY PRODUCTS.

THE SOFTWARE LICENSE AND LIMITED WARRANTY FOR THE ACCOMPANYING PRODUCT ARE SET FORTH IN THE INFORMATION PACKET THAT SHIPPED WITH THE PRODUCT AND ARE INCORPORATED HEREIN BY THIS REFERENCE. IF YOU ARE UNABLE TO LOCATE THE SOFTWARE LICENSE OR LIMITED WARRANTY, CONTACT YOUR CISCO REPRESENTATIVE FOR A COPY.

The following information is for FCC compliance of Class A devices: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio-frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference, in which case users will be required to correct the interference at their own expense.

Modifications to this product not authorized by Cisco could void the FCC approval and negate your authority to operate the product.

The Cisco implementation of TCP header compression is an adaptation of a program developed by the University of California, Berkeley (UCB) as part of UCB's public domain version of the UNIX operating system. All rights reserved. Copyright © 1981, Regents of the University of California.

NOTWITHSTANDING ANY OTHER WARRANTY HEREIN, ALL DOCUMENT FILES AND SOFTWARE OF THESE SUPPLIERS ARE PROVIDED "AS IS" WITH ALL FAULTS. CISCO AND THE ABOVE-NAMED SUPPLIERS DISCLAIM ALL WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING, WITHOUT LIMITATION, THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT OR ARISING FROM A COURSE OF DEALING, USAGE, OR TRADE PRACTICE.

IN NO EVENT SHALL CISCO OR ITS SUPPLIERS BE LIABLE FOR ANY INDIRECT, SPECIAL, CONSEQUENTIAL, OR INCIDENTAL DAMAGES, INCLUDING, WITHOUT LIMITATION, LOST PROFITS OR LOSS OR DAMAGE TO DATA ARISING OUT OF THE USE OR INABILITY TO USE THIS MANUAL, EVEN IF CISCO OR ITS SUPPLIERS HAVE BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

CCDE, CCENT, CCSI, Cisco Eos, Cisco HealthPresence, Cisco IronPort, the Cisco logo, Cisco Nurse Connect, Cisco Pulse, Cisco SensorBase, Cisco StackPower, Cisco StadiumVision, Cisco TelePresence, Cisco Unified Computing System, Cisco WebEx, DCE, Flip Channels, Flip for Good, Flip Mino, Flipshare (Design), Flip Ultra, Flip Video, Flip Video (Design), Instant Broadband, and Welcome to the Human Network are trademarks; Changing the Way We Work, Live, Play, and Learn, Cisco Capital, Cisco Capital (Design), Cisco:Financed (Stylized), Cisco Store, Flip Gift Card, and One Million Acts of Green are service marks; and Access Registrar, Aironet, AllTouch, AsyncOS, Bringing the Meeting To You, Catalyst, CCDA, CCDP, CCIE, CCIP, CCNA, CCNP, CCSP, CCVP, Cisco, the Cisco Certified Internetwork Expert logo, Cisco IOS, Cisco Lumin, Cisco Nexus, Cisco Press, Cisco Systems, Cisco Systems Capital, the Cisco Systems logo, Cisco Unity, Collaboration Without Limitation, Continuum, EtherFast, EtherSwitch, Event Center, Explorer, Follow Me Browsing, GainMaker, iLYNX, IOS, iPhone, IronPort, the IronPort logo, Laser Link, LightStream, Linksys, MeetingPlace, MeetingPlace Chime Sound, MGX, Networkers, Networking Academy, PCNow, PIX, PowerKEY, PowerPanels, PowerTV, PowerTV (Design), PowerVu, Prisma, ProConnect, ROSA, SenderBase, SMARTnet, Spectrum Expert, StackWise, WebEx, and the WebEx logo are registered trademarks of Cisco Systems, Inc. and/or its affiliates in the United States and certain other countries.

All other trademarks mentioned in this document or website are the property of their respective owners. The use of the word partner does not imply a partnership relationship between Cisco and any other company. (0910R)

Any Internet Protocol (IP) addresses and phone numbers used in this document are not intended to be actual addresses and phone numbers. Any examples, command display output, network topology diagrams, and other figures included in the document are shown for illustrative purposes only. Any use of actual IP addresses or phone numbers in illustrative content is unintentional and coincidental.

Cisco HealthPresence is intended for use by licensed health care professionals with those certain independent third party medical devices listed as compatible in the Cisco Health Presence Instructions for Use. The user is to refer to the third party instructions for use concerning any further information about the use of those certain medical devices. Cisco HealthPresence is not intended to perform real-time, active, or online patient monitoring, and does not transmit or display any real-time data that is intended to alert the user of alarms or other conditions that require immediate action or response. The user is advised that Cisco Systems makes no representations or warranties concerning jurisdictional requirements related to the licensed practice of medicine or health care using Cisco HealthPresence.

Cisco HealthPresence Solution Design Guide
© 2011 Cisco Systems, Inc. All rights reserved.

CHAPTER 1**Design Considerations for Systems Using Neurosynaptic Medical Devices 1-1**

Objectives and Audience 1-1

Cisco HealthPresence Intended Uses 1-1

Design Considerations When Using Neurosynaptic Medical Devices 1-2

Cisco HealthPresence-Supported Medical Devices and Peripherals 1-3

Neurosynaptic Devices 1-3

Bandwidth Requirements 1-4

Video Endpoint Bandwidth 1-4

Medical Device Bandwidth 1-4

Security 1-4



CHAPTER 1

Design Considerations for Systems Using Neurosynaptic Medical Devices

Revised: November 18, 2011

Objectives and Audience

This addendum describes the design considerations specific to an enterprise using the Neurosynaptic medical devices¹. This addendum must be used in conjunction with the Cisco HealthPresence™ Solution Design Guide, which provides an overview of the features and functions of Cisco HealthPresence and describes the options and best practices for designing a Cisco HealthPresence implementation, including:

- selecting a deployment model
- choosing the features you want to use
- designing the video conferencing portion of the solution
- determining bandwidth requirements and implementing quality of service
- determining security requirements and defining your security policy
- designing a Cisco HealthPresence room
- selecting Cisco service options

Almost all of the content of the Cisco HealthPresence Solution Design Guide applies to any implementation, regardless of the medical devices. This addendum covers the differences, which includes the devices supported, a minor security difference and a minor bandwidth difference. The target audience is Cisco HealthPresence solution planners and designers for a medical facility that plans to use Neurosynaptic medical devices.

Cisco HealthPresence Intended Uses

The Cisco HealthPresence Device consists of Attendant and Provider endpoints running the Cisco HealthPresence-Connect software.

1. Neurosynaptic devices are not approved for use in the United States and several other countries. For a list of countries where these devices are available for use, contact Neurosynaptic Communications Pvt. Ltd.

The Cisco HealthPresence Enterprise Server at a data center provides a communications link between HealthPresence endpoints over private or virtual private network connections having securement. Utilizing the network connections with securement, the Cisco HealthPresence-Connect software transmits data acquired from third-party medical devices at an Attendant endpoint to Provider endpoints located at remote facilities.

Compatible Third Party Devices include medical devices from AMD Telemedicine (covered in the Cisco HealthPresence Solution Design Guide) and the following medical devices from Neurosynaptic:

- ReMeDi ECG Probe (ReMeDi-ECGCAB-001)
- ReMeDi Temperature Probe (ReMeDi-TEMPCAB-001)
- ReMeDi Stethoscope Probe, (ReMeDi-STETHCAB-001)
- ReMeDi NiBP Probe, (ReMeDi-BPCAB-001)

The medical devices attach to the ReMeDi Medical Data Acquisition Unit (ReMeDi-MDAU), which in turn connects to the Cisco HealthPresence Attendant Appliance.

Additional network connections with securement, separate from those used to transmit the data from medical devices, provide video conferencing capabilities between an Attendant and a Provider endpoint.

The Cisco HealthPresence Connect software transmits data from a Cisco HealthPresence Attendant Appliance to a Cisco HealthPresence Provider Appliance, according to the intended uses of the third party medical devices.


Note

The Cisco HealthPresence-Connect software is not intended to perform real-time, active, or online patient monitoring, and does not transmit or display any real-time data that is intended to alert a physician of alarms or other conditions that require a physician's immediate action or response.

The user's infrastructure must meet the basic minimum specified requirements for Cisco HealthPresence to perform as intended. This includes, but is not limited to:

- adequate bandwidth
- appropriate video call control and routing equipment
- appropriate capacity routing and switching equipment
- proper software levels
- correct room design
- adequate physical and network security

Design Considerations When Using Neurosynaptic Medical Devices

The capabilities, components and architecture of the Cisco HealthPresence solution are the same regardless of which medical devices are attached to it; however, the attached medical devices determine the content that can be shared during an appointment and the tabs that appear in both the Provider and the Attendant screens. In the case of Neurosynaptic, that content is vitals data, ECG information and stethoscopic audio.

Cisco HealthPresence-Supported Medical Devices and Peripherals

Cisco HealthPresence is interoperable with both AMD Telemedicine² and Neurosynaptic devices³. Cisco HealthPresence also supports the Hi-Tech E-Pen (not a medical device) to enable electronic prescriptions.

The specific Neurosynaptic medical devices that are supported are described next.

Neurosynaptic Devices

In some countries outside of the United States, Cisco HealthPresence supports the Neurosynaptic Multi-parameter Data Acquisition Unit, referred to as NS MDAU. This device, shown in [Figure 1-1](#), provides the following measurements: 12-channel ECG, NiBP, Auscultation Sounds (Stethoscope), temperature and heart rate.

Figure 1-1 Neurosynaptic Multi-parameter Data Acquisition Unit



Bandwidth Requirements

Video Endpoint Bandwidth

The bandwidth for video endpoints is the same regardless of the medical devices attached to the Cisco HealthPresence Appliance. Refer to the Cisco HealthPresence Solution Design Guide for this information.

Medical Device Bandwidth

The Neurosynaptic medical devices require a minimum of 96Kbps and up to 191Kbps if all medical devices are transmitting concurrently. The breakdown of bandwidth is:

- ECG: 96 Kbps
- Stethoscope: 64 Kbps
- Temperature: 8 Kbps

2. AMD devices are not covered in this addendum.

3. Neurosynaptic devices are not supported in the United States.

- Blood Pressure: 8 Kbps
- Pulse Oximeter: 4.8 Kbps

Security

With the Neurosynaptic Device Aggregator the audio stream is transmitted via TCP and is not encrypted. This audio stream cannot be associated with a patient and all patient data is encrypted.